




provides an exceptional value for everyday oximetry needs. This compact, stand-alone monitor delivers fast, reliable SpO₂, pulse rate and pulse strength measurements complete with a plethysmographic waveform, and up to 30 hours of on-screen trending. Plus, the bright liquid crystal display, easy-to-use menu, and audible and visual alarms enhance clinical efficiency and quality of care for patients ranging from neonate to adult.

For more information visit smiths-medical.com/bci.

PATIENT MONITORING



Technical Specifications

DISPLAYS

- **Display**
Backlit LCD
- **Resolution**
240 x 64 pixels
- **Waveform**
Plethysmogram
- **Graphic Trends**
SpO₂, pulse rate (heart rate), pulse signal strength and respiration
- **SpO₂**
SpO₂ value, high and low alarm limits
- **Pulse Rate**
Pulse rate value, high and low alarm limits
- **Pulse Strength**
Logarithmically scaled 8-segment bar graph

INDICATORS

- **Low Battery**
Low battery icon when 10 minutes of battery life remain
- **Audio Status**
Yellow LED and displayed symbol indicates alarm audio paused (2 min) or indefinite
- **SpO₂ Sensor**
Message indicates sensor problem
- **SpO₂ Search**
Message indicates monitor is searching for a pulse
- **Small Pulse**
Message indicates oximeter signal strength is low
- **Charge**
Green CHARGE LED indicates external power and battery is charging
- **Lost Pulse**
Message indicates monitor has lost the signal

SPO₂

- **Range**
0 to 100% SpO₂ (functional)
- **Accuracy**
Adult: ± 2% at 70-100%
Neonate: ±3% at 70-100% SpO₂
<70% unspecified

- SpO₂, pulse rate and pulse strength measurements.
- Backlit liquid crystal display.
- Plethysmographic waveform.
- Up to 30 hours of on screen trended data for all parameters.
- Easy-to-operate menu structure.

Averaging

- 4, 8, or 16 pulse beat average
- **Alarm Range**
High: 50 to 100% and OFF (1% steps)
Low: 50 to 99% and OFF (1% steps)

PULSE RATE

- **Range**
30-254 bpm
- **Accuracy**
±2 bpm or ±2% (whichever is greater) at 30-254 bpm
- **Alarm Range**
High: 30-254 bpm and OFF (1 bpm steps)
Low: 30-254 bpm and OFF (1 bpm steps)
- **Averaging**
8 or 16 second average

POWER REQUIREMENTS

Two voltage ranges selected by external switch
90-125 VAC, 50/60 Hz
200-240 VAC, 50/60 Hz

BATTERY LIFE

Standard lead-acid 6 VDC battery:
Internal rechargeable
approximately 4 hrs of continuous use (new)
approximately 7 hrs of charge time

PRINTER OUTPUT

Serial infrared link—real time patient data RS232C trend table output

DIMENSIONS

8" W x 6.5" H x 5" D
(20.3 cm W x 16.5 cm H x 12.7 cm D)

WEIGHT

4.5 lbs. (2kg)

ENVIRONMENTAL

- **Operating Temperature**
32°F to 131°F (0°C to 55°C)
- **Storage Temperature**
-40°F to 167°F (-40° to 75°C)
- **Relative Humidity** (non-condensing)
10-95% (storage)
15-95% (operating)

- Audible and visual alarms.
- AC or battery power.
- Monitor patients ranging from neonate to adult.
- Genuine BCI® reusable and disposable sensors available.

STANDARDS

EN60601-1 (IEC 601-1)
UL 2601-1 and CAN/CSA 601.1
CAN/CSA NRTL/C

REGULATORY COMPLIANCE

ISO 9001/ISO 13485

CATALOG NUMBER

- **3180**
Pulse Oximeter

Information subject to change without notice. Please read the Instructions for Use supplied with the product for detailed instructions, warnings and cautions.

Q-Line for U.S.A. customers: call 1-800-348-1842 for clinical application information.

For more information, please call Smiths Medical PM, Inc., at 262-542-3100 or 800-558-2345, or your Smiths Medical distributor.

Smiths Medical PM, Inc.

Patient Monitoring and Ventilation

N7W22025 Johnson Drive, Waukesha, Wisconsin 53186 USA

Phone: 262-542-3100 Fax: 262-542-0718 Toll-Free: 800-558-2345

www.smiths-medical.com/bci

Smiths Medical, part of the global technology business Smiths Group

BCI and the Smiths Medical design mark are trademarks of the Smiths Medical family of companies.

The symbol ® indicates the trademark is registered in the U.S. Patent and Trademark Office and certain other countries.

All other names and marks mentioned are the trade names, trademarks or service marks of their respective owners.

©2008 Smiths Medical family of companies. All rights reserved. BCI 3180 - PM195607EN Rev.02 06/08

smiths medical